

Codes of Teleopto Cannules

Code in product name	Meaning	Select Your need			Notes
c	color of LED	B: 470nm(Blue)	G: 525nm(Green)	Y: 590nm(Yellow)	Please aske SANSUIDO if need other wavelength
l	Length of cannule	Please specify length you need (in mm unit)			
d	diameter of optical fiber	250: (φ 250nm)	500: (φ 500nm)	750: (φ 750nm)	
i	interval of between two cannules	Please specify length you need (in mm unit)			

Reference Datas for Your Select --- Optoical Power out-put

Led color / Cannule diameter	out-put	Notes
LED B: 470nm / Cannule φ250nm	4.5mW (≒99.5mW/mm ²)	Using Teler-2-P + TeleLC
LED B: 470nm / Cannule φ500nm	13.0mW (≒71.8mW/mm ²)	
LED G: 525nm / Cannule φ250nm	1.4mW (≒31.1mW/□)	
LED G: 525nm / Cannule φ500nm	3.0mW (≒16.6mW/□)	
LED Y: 590nm / Cannule φ250nm	1.0mW (≒22.2mW/□)	
LED Y: 590nm / Cannule φ500nm	2.8mW (≒15.5mW/□)	